



ACNS Recommendations Equipment Maintenance and Cleaning

NOTE: The situation surrounding COVID-19 is evolving rapidly and there are unique aspects at the institutional and local levels. Guidance provided by your institution, state and federal authorities should always take precedence over the recommendations provided here.

Equipment Maintenance and Cleaning

1. Use antiseptic wipe to clean all surfaces of the NDT equipment that has entered any COVID+/PUI patient room
2. Consider clear plastic bags to cover EEG equipment in COVID+/PUI rooms
3. [See University of Maryland equipment cleaning protocol](#)